


CRE-2110 POWER OPTIONS TECH DRAWING




POWER SWIVEL SEAT ONLY
(USES MIDDLE MOUNTING HOLES,
LOWEST USABLE FOOTREST POSITION ON CRE-2110)
(SEE DETAIL A)




ILS-00937 SHEET 1 OF 1 REV. 1 (7530)(09/16/22)(BHM)



POWER FOLDING FOOTREST ONLY
(USES TOP MOUNTING HOLES,
LOWEST USABLE FOOTREST POSITION ON CRE-2110)
(SEE DETAIL B)



POWER SWIVEL SEAT AND POWER FOLDING FOOTREST (USES TOP MOUNTING HOLES, LOWEST USABLE FOOTREST POSITION ON CRE-2110) (SEE DETAIL B)



NOTE: WHEN THE BOTTOM FOOTREST MOUNTING HOLES ARE USED, THE SEAT WILL BE IN ITS HIGHEST POSITION.



©Bruno Independent Living Aids, Inc.